

# Non-Invasive Cholesterol Measurement System CLE-100



■ Sample measurement section (ATR prism)



■ Main instrument

The CLE-100 makes a simple non-invasive measurement of the cholesterol ester density in a subject's skin and reports the total cholesterol and LDL values in the blood. To make a measurement; place the patients arm on the measurement area for approximately one minute\*. The results are quickly calculated and reported.

\*measurement time may vary according to conditions.

## [ Features ]

- No blood sample required, the cholesterol ester quantity is measured directly.
- No reagent or pretreatment required.
- Compact, desk-top instrument
- Quick measurement using a combination of advanced infra-red detection and a high accuracy band pass filter.
- Measurement and data processing are controlled by PC

## [ Use ]

- Measurement of cholesterol ester and lipids in the skin



(Note) Patent applied

## CLE-100 - Non-Invasive Cholesterol Measurement System

[Example of measurement]



■ Measurement in progress (example)



■ FT-IR spectrum of subject measurement

**Red:** Spectrum for an individual with normal cholesterol levels

**Black:** Spectrum for an individual with high cholesterol levels

Cholesterol ester peak located at 1742cm-1.

■ Data shows high correlation with LDL levels in the blood

Absorption value at 1742cm-1, R=0.984

Absorption value at 1170cm-1, R=0.985

Specifications	CLE-100	※Product specification are subject to change without notice
Light source	Filament type IR source	
Measurement	ATR prism (ZnSe)	
Spectral system	High accuracy band pass filter (automatic filter change)	
Detector	High sensitivity DTGS	
PC for control and data processing	PC	
Software	OS: Windows 7 or Windows XP	
Dimensions	Main unit: 350(W)×250(D)×270(H), Accessories 320(W)×280(D)×250(H)	
Power source and power consumption	AC100V 50/60Hz · 300W (Max)	

This system was developed with a "2009 product development support subsidy for small and medium-sized businesses" (which supports businesses such as trial manufacture development)

Distributor

Manufacturer  
**Photoscience Incorporated**  
 1-D Kimura-bldg 492-1 Katakura  
 Hachioji, Tokyo 192-0914, Japan  
 Tel : +81-42-649-1447 Fax : +81-42-649-1455  
 URL <http://photoscience.co.jp>