

Photoscience

# Micro Raman Spectrometer EGR-100/300



The EGR-100 Micro Raman system offers an excellent cost performance ratio and despite its compact size (only 60cm width x 40 cm depth and only 25kg\*), the EGR-100 delivers performance equal to much larger Raman spectrometers.


The EGR-100 is fitted with an exceptional anti-vibration mechanism, the detector is a thermoelectric cooled CCD (down to -70 deg C). The advanced compact 532nm green laser is built into the optical bench and a mechanical interlock guarantees user safety\*\*.

\*Outline size

\*\*This instrument is in "class 1 laser" category. JIS C 6802: 2005

### [ Features ]

- Compact, light weight and low price but with exceptional performance
- Anti vibration microscopic system
- Compact spectrograph with advanced thermoelectric cooled CCD detector and high throughput optical system
- Advanced compact green laser (mechanical safety interlock)
- Optional horizontal irradiation of laser light and polarization measurement



Anti-vibration system

### [ Use ]

- Analysis of micro samples
- Forensic analysis of micro foreign material
- Structural analysis of DLC
- Structural analyses of polymer
- Liquid analysis
- Structural analysis using polarimetry (option)

**コメントの追加 [ 1]:** monochromator  
光学系をここでは、spectroscope としておりますが、monochromator や spectrograph など、適切な言葉に置き換えてください。

**コメントの追加 [ 2]:** check this translation for technical correctness この翻訳ですが、水平照射に vertical を訳者が使いましたが、horizontal が適切でしょうか？

**コメントの追加 [ 3]:** High molecular weight or complex structure, not sure which?  
ここでは、訳者より上記にありますように、どちらが適切がわかりませんでしたので、適切な方を選択して下さい。

Photoscience


## Photoscience Micro Raman Spectrometer EGR-100

### [Example of Raman measurement]

<Measurement conditions> Laser 532nm, output 9mW (sample position) and object lens  $\times 100$  (NA/0.8), 100  $\mu\text{m}$  slit width (band width approximately 0.7nm: approximately  $24\text{cm}^{-1}$ )


#### ★ Artificial diamond crystal

Green spot in the center is the laser light




For observation of the entire sample, a low magnification ratio lens  $\times 5$  is used.

[Measurement display of Andor CCD]




1 second measurement time  
[Raman spectrum of diamond]



Strong peak at  $1332\text{cm}^{-1}$  is observed.


#### ★ Polystyrene board

Measurement time: 30 secs




#### ★ Carbon ('HB' pencil lead with combination of nano-diamond)

Measurement time: 30 seconds



In addition to the D and G bands, secondary Raman scattering peaks DR2 and G' are also observed on the high wavenumber side.



■ Configuration

### Specifications

### EGR-100

※Product specification are subject to change without notice (2010. 4)

Excitation Laser	532nm green laser 55mW type (various laser enabled)
Sample chamber and sample stage	Sample stage: XYZ micro-pitch mechanism by micrometer
Objective lens system	Objective lens: $\times 100$ , $\times 50$ , $\times 20$ , $\times 5$ turret fitted as standard
Sample observation	Sample observation and the laser irradiation position verification by CCD camera
Laser irradiation system	Vertical irradiation (backscatter optical measurement) Horizontal irradiation ( $90^\circ$ scattered light measurement: Option)
Spectrograph	ANDOR corporation, Shamrock SR16 (Focal length 163mm and F/3.6) (standard)
Detector	Standard: ANDOR corporation DV401-F1 type CCD (Various CCD options)
PC for control and data processing	Laptop PC; OS: WindowsXP® or Vista® (2 or more USB terminals required)
Software	CCD control software and spectrum data processing software
Instrument dimensions and weight	Approx. 600(W) $\times$ 400(D) $\times$ 470(H) excluding projections, Approx. 25kg (Spectroscope approx. 3.5kg and CCD approx. 2.5kg)
Power source and power consumption	AC100V 50/60Hz ; Approx. 270W (max)

#### Manufacturer

#### Photoscience Incorporated

1-D Kimura-bldg 492-1 Katakura

Hachioji, Tokyo 192-0914, Japan

Tel : +81-42-649-1447 Fax : +81-42-649-1455

URL <http://photoscience.co.jp>